



Material Data Sheet

General Information

Headquarters

TELIX Pharmaceuticals Ltd
55 Flemington Road, Level 4
North Melbourne, VIC 3051
Australia

Manufacturing Facility

TELIX Innovations SA
Rue de Hermée, 255
4040 Herstal
Belgium

Material Ordering Description and Overview

For EZAG Generator

Optimized for 5 mL
eluate of 0.1 M HCl



Box of 5 Sterile Acetate Buffer Glass Vials, Config. A

ITEM CODE
900018

PRODUCT REFERENCE
BOX-GSV-BF-AC-A

SHELF LIFE
36 months

For IRE Generator

Optimized for 1.1 mL
eluate of 0.1 M HCl



Box of 5 Sterile Acetate Buffer Glass Vials, Config. B

ITEM CODE
900019

PRODUCT REFERENCE
BOX-GSV-BF-AC-B

SHELF LIFE
36 months

Sterile Vacuumed Glass Vial



Box of 5 Sterile Vacuumed Glass Vials Compatible with both generators

ITEM CODE
900017

PRODUCT REFERENCE
BOX-SVV-10G

SHELF LIFE
24 months

Storage conditions for all materials above: Do not freeze.